

Product Datasheet

Beta galactosidase protein (orb371914)

Description	Recombinant E. coli Beta-galactosidase protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Application notes	(M443L, C500S)
Protein Sequence	Met1-Lys1024(12-41AA deletion)
Source	E.coli
Uniprot ID	P00722
Expiration Date	6 months from date of receipt.