

Product Datasheet

Human Left-right Determination Factor 2 protein (orb371920)

Description	Recombinant human Left-right Determination Factor 2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	(N-6His)
Protein Sequence	Phe78-Pro366
Source	Human Cells
Uniprot ID	O00292
Expiration Date	6 months from date of receipt.