

Product Datasheet

TGF beta Receptor 2 antibody (orb371959)

Description

TGF beta Receptor 2 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human TGF beta Receptor 2. Please contact us for the exact immunogen sequence. The peptide is available as orb374981.

Target

TGF beta Receptor 2

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

56 kDa

Uniprot ID

[P38438](#), [P37173](#), [Q62312](#)

Entrez



[7048](#)

Dilution Range


IHC-P: 1:100-500, WB: 1:100-1000, IF/ICC: 1:100-500 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



IHC-P staining of rat pancreas tissue using TGF beta Receptor 2 antibody.



IF analysis of rat brain tissue using TGF beta Receptor 2 antibody.



IHC-P image of rat brain tissue using TGF beta Receptor 2 antibody.