

Product Datasheet

TNF alpha antibody (orb371962)

Description

TNF alpha antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human TNF alpha. Please contact us for the exact immunogen sequence. The peptide is available as orb374984.

Target

TNF alpha

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

27 kDa

Uniprot ID

[P06804](#), [P16599](#), [P01375](#)

NCBI

[000594](#), [000585](#)

Entrez



[7124](#)

Dilution Range


IF/ICC: 1:100-800, WB: 1:100-2000, IHC-P: 1:100-800 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



IHC-P image of human breast tumour using...



Immunofluorescence analysis of rat skin ...



Immunohistochemical staining of rat skin...