

Product Datasheet

Granzyme B antibody (orb372089)

Description	Rabbit polyclonal antibody to Granzyme B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P: 1:10 - 1:50
Expiration Date	12 months from date of receipt.