

Product Datasheet

NLRP1 antibody (FITC) (orb37220)

Description	NLRP1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 1 protein (1-146AA)
Target	NLRP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9C000
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.