

Product Datasheet

Pro-GRP antibody (Coating) (orb372234)

Description	Mouse monoclonal antibody to Pro-GRP (Coating)
Species/Host	Mouse
Conjugation	Unconjugated
Target	Pro-GRP
Preservatives	PBS pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Application notes	Suggested applications: EIA, CLIA, RIA, IHC, WB
Isotype	IgG
Clonality	Monoclonal
Source	Ascites
Expiration Date	12 months from date of receipt.