

## Product Datasheet

### ICA69 antibody (orb372267)

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Description</b>       | Mouse monoclonal antibody to ICA69         |
| <b>Species/Host</b>      | Rabbit                                     |
| <b>Conjugation</b>       | Unconjugated                               |
| <b>Target</b>            | ICA69                                      |
| <b>Note</b>              | For research use only                      |
| <b>Application notes</b> | Suggested applications: LFA, EIA, CLIA, WB |
| <b>Isotype</b>           | IgG                                        |
| <b>Clonality</b>         | Polyclonal                                 |
| <b>Expiration Date</b>   | 12 months from date of receipt.            |