



Product Datasheet

PEDF antibody (orb372324)


Description	PEDF antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant protein of human PEDF
Target	PEDF
Form/Appearance	PBS with 0.05% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified by protein G.
MW	47 kDa
Uniprot ID	P36955, P97298
NCBI	002606
Entrez	5176
Dilution Range	WB:1:1000-1:2000, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



IHC-P image of mouse skin tissue using a...



Western blot analysis of Pedf recombinan...



Immunohistochemical staining of rat lym...