



Product Datasheet

CD133 antibody (orb372328)


Description	CD133 antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein of human CD133
Target	CD133
Form/Appearance	PBS with 0.05% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified by protein G.
MW	30.8 kDa
Uniprot ID	Q8CJ52 , O54990 , O43490
NCBI	001139319 , 001145847
Entrez	8842
Dilution Range	WB: 1:200-2000, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



WB analysis of human lung cancer (lane 1...



Immunohistochemical staining of mouse th...



Western blot analysis of CD133-E1 recomb...