



Product Datasheet

Dnmt3b antibody (orb372330)


Description	Dnmt3b antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Dnmt3b. Please contact us for the exact immunogen sequence. The peptide is available as orb374989.
Target	Dnmt3b
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	96 kDa
Uniprot ID	Q1LZ51, Q9UBC3, O88509
NCBI	001207055, 001193984
Entrez	1789
Dilution Range	WB: 1:100-500, IHC-P: 1:50-400
Expiration Date	12 months from date of receipt.



IHC-P
staining of
rat lung
tissue using
...



IHC-P
staining of
rat lung
tissue using
...



Western blot
analysis of
rat heart
(lane...