

Product Datasheet

HBV polymerase (455-463) (orb372349)

Description	Peptide to HBV polymerase (455-463)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Leu-Ser-Arg-Tyr-Val-Ala-Arg- Leu
Formula	$C_{46}H_{79}N_{15}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1034.24
Expiration Date	12 months from date of receipt.