

Product Datasheet

HBV polymerase (502-510) (orb372350)

Description	Peptide to HBV polymerase (502-510)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Leu-His-Leu-Tyr-Ser-His-Pro-Ile
Formula	$C_{53}H_{82}N_{14}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1107.33
Expiration Date	12 months from date of receipt.