

Product Datasheet

Helodermin (orb372364)

Description	Peptide to Helodermin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Ser-Asp-Ala-Ile-Phe-Thr-Glu-Glu-Tyr-Ser-Lys-Leu-Leu-Ala-Lys-Leu-Ala-Leu-Gln-Lys-Tyr-Leu-Ala-Ser-Ile-Leu-Gly-Ser-Arg-Thr-Ser-Pro-Pro-Pro-NH2
Formula	C ₁₇₆ H ₂₈₃ N ₄₅ O ₅₁
Purity	> 95%
Source	Synthetic
MW	3845.49
CAS Number	89468-62-2
Expiration Date	12 months from date of receipt.