

Product Datasheet

Hemokinin 1 (orb372368)

Description	Peptide to Hemokinin 1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Gly-Lys-Ala-Ser-Gln-Phe-Phe-Gly-Leu-Met-NH ₂
Formula	C ₅₄ H ₈₄ N ₁₄ O ₁₄ S
Purity	> 95%
Source	Synthetic
MW	1185.42
CAS Number	491851-53-7
Expiration Date	12 months from date of receipt.

