

Product Datasheet

Hepatitis Virus C NS3 Protease Inhibitor 1 (orb372371)

Description	Peptide to Hepatitis Virus C NS3 Protease Inhibitor 1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ac-Asp-Glu-Met-Glu-Glu-Cys
Formula	$C_{29}H_{45}N_6O_{16}S_2$
Purity	> 95%
Source	Synthetic
MW	796.80
Expiration Date	12 months from date of receipt.