

Product Datasheet

Histone H3 (116-136) (orb372380)

Description	Peptide to Histone H3 (116-136)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Arg-Val-Thr-Ile-Met-Pro-Lys-Asp-Ile-Gln-Leu-Ala-Arg-Arg-Ile-Arg-Gly-Glu-Arg-Ala
Formula	$C_{107}H_{195}N_{39}O_{28}S$
Purity	> 95%
Source	Synthetic
MW	2508.06
Expiration Date	12 months from date of receipt.