

Product Datasheet

Histone H3 (1-25) amide (orb372382)

Description	Peptide to Histone H3 (1-25) amide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ala-Arg-Thr-Lys-Gln-Thr-Ala-Arg-Lys-Ser-Thr-Gly-Gly-Lys-Ala-Pro-Arg-Lys-Gln-Leu-Ala-Thr-Lys-Ala-Ala-NH ₂
Formula	C ₁₁₀ H ₂₀₂ N ₄₂ O ₃₂
Purity	> 95%
Source	Synthetic
MW	2625.10
Expiration Date	12 months from date of receipt.