

Product Datasheet

HIV-1 gag Protein p24 (137-154) (orb372388)

Description	Peptide to HIV-1 gag Protein p24 (137-154)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Leu-Asn-Lys-Ile-Val-Arg-Met-Tyr-Ser-Pro-Thr-Ser-Ile-Leu-Asp-Ile-Arg
Formula	$C_{92}H_{158}N_{26}O_{26}S$
Purity	> 95%
Source	Synthetic
MW	2076.51
Expiration Date	12 months from date of receipt.

