

Product Datasheet

HIV-1 tat Protein (49-57) (orb372394)

Description	Peptide to HIV-1 tat Protein (49-57)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Lys-Lys-Arg-Arg-Gln-Arg-Arg-Arg
Formula	$C_{53}H_{106}N_{30}O_{11}$
Purity	> 95%
Source	Synthetic
MW	1339.63
CAS Number	123251-89-8
Expiration Date	12 months from date of receipt.

