

Product Datasheet

HIV-gp120-Antigenic Peptide (orb372399)

Description	Peptide to HIV-gp120-Antigenic Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Cys-Gly-Lys-Ile-Glu-Pro-Leu-Gly-Val-Ala-Pro-Thr-Lys-Ala-Lys-Arg-Arg-Val-Val-Gln-Arg-Glu-Lys-Arg
Formula	$C_{117}H_{211}N_{41}O_{31}S$
Purity	> 95%
Source	Synthetic
MW	2720.30
CAS Number	198636-94-1
Expiration Date	12 months from date of receipt.

