

## Product Datasheet

### HJ Inhibitor Peptide 2 (orb372406)

|                         |                                              |
|-------------------------|----------------------------------------------|
| <b>Description</b>      | Peptide to HJ Inhibitor Peptide 2            |
| <b>Form/Appearance</b>  | Lyophilized powder                           |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year. |
| <b>Note</b>             | For research use only                        |
| <b>Protein Sequence</b> | Lys-Trp-Trp-Cys-Arg-Trp                      |
| <b>Formula</b>          | $C_{48}H_{61}N_{13}O_7S$                     |
| <b>Purity</b>           | > 95%                                        |
| <b>Source</b>           | Synthetic                                    |
| <b>MW</b>               | 964.16                                       |
| <b>Expiration Date</b>  | 12 months from date of receipt.              |