

Product Datasheet

Human CMV Assemblin Protease Substrate (orb372416)

Description	Peptide to Human CMV Assemblin Protease Substrate
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Arg-Gly-Val-Val-Asn-Ala-Ser-Ser-Arg- Leu-Ala-Lys
Formula	$C_{52}H_{96}N_{20}O_{16}$
Purity	> 95%
Source	Synthetic
MW	1257.47
Expiration Date	12 months from date of receipt.

