

Product Datasheet

Human IgE Pentapeptide HEPP (orb372418)

Description	Peptide to Human IgE Pentapeptide HEPP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Ser-Asp-Pro-Arg
Formula	$C_{22}H_{36}N_8O_{11}$
Purity	> 95%
Source	Synthetic
MW	588.58
CAS Number	62087-72-3
Expiration Date	12 months from date of receipt.

