

Product Datasheet

Hydrin 2 (orb372422)

Description	Peptide to Hydrin 2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Cys-Tyr-Ile-Gln-Asn-Cys-Pro-Arg-Gly-Gly (Disulfide bridge Cys1-Cys6)
Formula	$C_{45}H_{69}N_{15}O_{14}S_2$
Purity	> 95%
Source	Synthetic
MW	1108.27
CAS Number	122842-55-1
Expiration Date	12 months from date of receipt.

