

Product Datasheet

Hypertrehalosaemic Neuropeptide (orb372426)

Description	Peptide to Hypertrehalosaemic Neuropeptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Pyr-Leu-Thr-Phe-Ser-Ser-Gly-Trp-Gly-Asn-NH ₂
Formula	C ₄₉ H ₆₇ N ₁₃ O ₁₅
Purity	> 95%
Source	Synthetic
MW	1078.17
Expiration Date	12 months from date of receipt.

