

Product Datasheet

Hypertrehalosaemic Neuropeptide (orb372427)

Description	Peptide to Hypertrehalosaemic Neuropeptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Pyr-Val-Asn-Phe-Ser-Pro-Gly-Trp-Gly-Thr-NH ₂
Formula	C ₅₀ H ₆₇ N ₁₃ O ₁₄
Purity	> 95%
Source	Synthetic
MW	1074.19
CAS Number	106018-36-4
Expiration Date	12 months from date of receipt.

