

Product Datasheet

IELLQAR Selectin Binding Peptide (orb372429)

Description	Peptide to IELLQAR Selectin Binding Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Glu-Leu-Leu-Gln-Ala-Arg
Formula	$C_{37}H_{67}N_{11}O_{11}$
Purity	> 95%
Source	Synthetic
MW	842.01
Expiration Date	12 months from date of receipt.

