

Product Datasheet

IGF-I (30-41) (orb372432)

Description	Peptide to IGF-I (30-41)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Gly-Tyr-Gly-Ser-Ser-Ser-Arg-Arg-Ala-Pro-Gln-Thr
Formula	$C_{51}H_{83}N_{19}O_{19}$
Purity	> 95%
Source	Synthetic
MW	1266.35
CAS Number	82177-09-1
Expiration Date	12 months from date of receipt.

