

Product Datasheet

IGF-II (33-40) (orb372434)

Description	Peptide to IGF-II (33-40)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ser-Arg-Val-Ser-Arg-Arg-Ser-Arg
Formula	$C_{38}H_{74}N_{20}O_{12}$
Purity	> 95%
Source	Synthetic
MW	1003.14
CAS Number	103244-41-3
Expiration Date	12 months from date of receipt.

