

## Product Datasheet

### IL-1a (223-250) (orb372436)

|                         |                                                                                                                 |
|-------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Peptide to IL-1a (223-250)                                                                                      |
| <b>Form/Appearance</b>  | Lyophilized powder                                                                                              |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year.                                                                    |
| <b>Note</b>             | For research use only                                                                                           |
| <b>Protein Sequence</b> | Phe-Phe-Tyr-Glu-Thr-His-Gly-Thr-Lys-Asn-Tyr-Phe-Thr-Ser-Val-Ala-His-Pro-Asn-Leu-Phe-Ile-Ala-Thr-Lys-Gln-Asp-Tyr |
| <b>Formula</b>          | $C_{158}H_{219}N_{37}O_{44}$                                                                                    |
| <b>Purity</b>           | > 95%                                                                                                           |
| <b>Source</b>           | Synthetic                                                                                                       |
| <b>MW</b>               | 3340.73                                                                                                         |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                                                                 |

