

Product Datasheet

IL-1b (208-240) (orb372438)

Description	Peptide to IL-1b (208-240)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Lys-Lys-Lys-Met-Glu-Lys-Arg-Phe-Val-Phe-Asn-Lys-Ile-Glu-Ile-Asn-Asn-Lys-Leu-Glu-Phe-Glu-Ser-Ala-Gln-Phe-Pro-Asn-Trp-Tyr-Ile-Ser-Thr
Formula	$C_{191}H_{292}N_{48}O_{51}S$
Purity	> 95%
Source	Synthetic
MW	4108.81
Expiration Date	12 months from date of receipt.