

Product Datasheet

IL-6 (88-121) (orb372439)

Description	Peptide to IL-6 (88-121)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Ile-Thr-Gly-Leu-Leu-Glu-Phe-Glu-Val-Tyr-Leu-Glu-Tyr-Leu-Gln-Asn-Arg-Phe-Glu-Ser-Ser-Glu-Glu-Gln-Ala-Arg-Ala-Val-Gln-Met-Ser-Thr-Lys
Formula	$C_{179}H_{281}N_{45}O_{58}S$
Purity	> 95%
Source	Synthetic
MW	4023.57
CAS Number	145990-81-4
Expiration Date	12 months from date of receipt.