

Product Datasheet

Insulin Receptor (1142-1153), pTyr1151 (orb372462)

Description	Peptide to Insulin Receptor (1142-1153), pTyr1151
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Arg-Asp-Ile-Tyr-Glu-Thr-Asp-Tyr-pTyr-Arg-Lys
Formula	$C_{72}H_{108}N_{19}O_{27}$
Purity	> 95%
Source	Synthetic
MW	1702.77
Expiration Date	12 months from date of receipt.