

Product Datasheet

Interleukin II Beta (44-56) (orb372467)

Description	Peptide to Interleukin II Beta (44-56)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Leu-Asn-Gly-Ile-Asn-Asn-Tyr-Lys-Asn-Pro-Lys-Leu
Formula	$C_{68}H_{113}N_{19}O_{19}$
Purity	> 95%
Source	Synthetic
MW	1500.77
Expiration Date	12 months from date of receipt.

