

Product Datasheet

KGF Receptor Peptide (orb372478)

Description	Peptide to KGF Receptor Peptide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	His-Ser-Gly-Ile-Asn-Ser-Ser-Asn-Ala-Glu-Val-Leu-Ala-Leu-Phe-Asn-Val-Thr-Glu-Met-Asp-Ala-Gly-Glu-Tyr
Formula	$C_{114}H_{174}N_{30}O_{42}S$
Purity	> 95%
Source	Synthetic
MW	2668.90
Expiration Date	12 months from date of receipt.

