

Product Datasheet

Kisspeptin-54 (27-54) (orb372483)

Description	Peptide to Kisspeptin-54 (27-54)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Pro-Ala-Pro-Gln-Gly-Ala-Val-Leu-Val-Gln-Arg-Glu-Lys-Asp-Leu-Pro-Asn-Tyr-Asn-Trp-Asn-Ser-Phe-Gly-Leu-Arg-Phe-NH ₂
Formula	C ₁₄₉ H ₂₂₆ N ₄₂ O ₃₉
Purity	> 95%
Source	Synthetic
MW	3229.72
Expiration Date	12 months from date of receipt.