

## Product Datasheet

### KK4A (orb372484)

|                         |                                                                 |
|-------------------------|-----------------------------------------------------------------|
| <b>Description</b>      | Peptide to KK4A                                                 |
| <b>Form/Appearance</b>  | Lyophilized powder                                              |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year.                    |
| <b>Note</b>             | For research use only                                           |
| <b>Protein Sequence</b> | Lys-Lys-Gly-Ser-Val-Val-Ile-Val-Gly-Arg-Ile-Val-Leu-Ser-Gly-Lys |
| <b>Formula</b>          | $C_{74}H_{138}N_{22}O_{19}$                                     |
| <b>Purity</b>           | > 95%                                                           |
| <b>Source</b>           | Synthetic                                                       |
| <b>MW</b>               | 1640.06                                                         |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                 |

