

Product Datasheet

Laminin (925-933) (orb372485)

Description	Peptide to Laminin (925-933)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Cys-Asp-Pro-Gly-Tyr-Ile-Gly-Ser-Arg
Formula	$C_{40}H_{62}N_{12}O_{14}S$
Purity	> 95%
Source	Synthetic
MW	967.10
CAS Number	110590-60-8
Expiration Date	12 months from date of receipt.

