

Product Datasheet

Laminin (929-933) (orb372486)

Description	Peptide to Laminin (929-933)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Tyr-Ile-Gly-Ser-Arg
Formula	$C_{26}H_{42}N_8O_8$
Purity	> 95%
Source	Synthetic
MW	594.70
CAS Number	110590-64-2
Expiration Date	12 months from date of receipt.

