

Product Datasheet

Leptin (93-105) (orb372492)

Description	Peptide to Leptin (93-105)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asn-Val-Ile-Gln-Ile-Ser-Asn-Asp-Leu-Glu-Asn-Leu-Arg
Formula	$C_{64}H_{110}N_{20}O_{23}$
Purity	> 95%
Source	Synthetic
MW	1527.70
CAS Number	200436-43-7
Expiration Date	12 months from date of receipt.

