

## Product Datasheet

### Leucokinin I (orb372494)

|                         |                                                                 |
|-------------------------|-----------------------------------------------------------------|
| <b>Description</b>      | Peptide to Leucokinin I                                         |
| <b>Form/Appearance</b>  | Lyophilized powder                                              |
| <b>Storage</b>          | Shipped at 4°C. Store at -20°C for one year.                    |
| <b>Note</b>             | For research use only                                           |
| <b>Protein Sequence</b> | Asp-Pro-Ala-Phe-Asn-Ser-Trp-Gly-NH <sub>2</sub>                 |
| <b>Formula</b>          | C <sub>41</sub> H <sub>52</sub> N <sub>11</sub> O <sub>12</sub> |
| <b>Purity</b>           | > 95%                                                           |
| <b>Source</b>           | Synthetic                                                       |
| <b>MW</b>               | 891.93                                                          |
| <b>CAS Number</b>       | 104600-89-7                                                     |
| <b>Expiration Date</b>  | 12 months from date of receipt.                                 |

