

Product Datasheet

Leucokinin II (orb372495)

Description	Peptide to Leucokinin II
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Pro-Gly-Phe-Ser-Ser-Trp-Gly-NH ₂
Formula	C ₃₉ H ₅₀ N ₁₀ O ₁₂
Purity	> 95%
Source	Synthetic
MW	850.90
CAS Number	104821-31-0
Expiration Date	12 months from date of receipt.

