

Product Datasheet

Leucokinin III (orb372496)

Description	Peptide to Leucokinin III
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Asp-Gln-Gly-Phe-Asn-Ser-Trp-Gly-NH ₂
Formula	C ₄₀ H ₅₂ N ₁₂ O ₁₃
Purity	> 95%
Source	Synthetic
MW	908.90
CAS Number	104958-71-6
Expiration Date	12 months from date of receipt.

