

Product Datasheet

Leucopyrokinin (orb372500)

Description	Peptide to Leucopyrokinin
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	PYR-Thr-Ser-Phe-Thr-Pro-Arg-Leu-NH ₂
Formula	C ₄₂ H ₆₆ N ₁₂ O ₁₂
Purity	> 95%
Source	Synthetic
MW	931.06
CAS Number	104052-00-8
Expiration Date	12 months from date of receipt.

