

Product Datasheet

Leucopyrokinin (4-8) (orb372501)

Description	Peptide to Leucopyrokinin (4-8)
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Phe-Thr-Pro-Arg-Leu-NH ₂
Formula	C ₃₀ H ₄₉ N ₉ O ₆
Purity	> 95%
Source	Synthetic
MW	631.78
CAS Number	104191-78-8
Expiration Date	12 months from date of receipt.

