

Product Datasheet

Leu-Valorphin-Arg (orb372502)

Description	Peptide to Leu-Valorphin-Arg
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Val-Val-Tyr-Pro-Trp-Thr-Gln-Arg
Formula	$C_{56}H_{84}N_{14}O_{13}$
Purity	> 95%
Source	Synthetic
MW	1161.38
CAS Number	135861-77-7
Expiration Date	12 months from date of receipt.

