

Product Datasheet

LL-37 Pentamide (orb372508)

Description	Peptide to LL-37 Pentamide
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Leu-Leu-Gly-Asn-Phe-Phe-Arg-Lys-Ser-Lys-Gln-Lys-Ile-Gly-Lys-Gln-Phe-Lys-Arg-Ile-Val-Gln-Arg-Ile-Lys-Asn-Phe-Phe-Arg-Asn-Leu-Val-Pro-Arg-Thr-Gln-Ser
Formula	$C_{208}H_{343}N_{65}O_{48}$
Purity	> 95%
Source	Synthetic
MW	4522.46
Expiration Date	12 months from date of receipt.