

Product Datasheet

LMP1 (156-164), IAL (orb372510)

Description	Peptide to LMP1 (156-164), IAL
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Ile-Ala-Leu-Tyr-Leu-Gln-Gln-Asn-Trp
Formula	$C_{55}H_{81}N_{13}O_{14}$
Purity	> 95%
Source	Synthetic
MW	1148.34
Expiration Date	12 months from date of receipt.