

Product Datasheet

LMP2 (419-427), TYG (orb372513)

Description	Peptide to LMP2 (419-427), TYG
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Protein Sequence	Thr-Tyr-Gly-Pro-Val-Phe-Met-Cys-Leu
Formula	$C_{48}H_{71}N_9O_{12}S_2$
Purity	> 95%
Source	Synthetic
MW	1030.27
Expiration Date	12 months from date of receipt.